Disposal and Recycling:

The Evans Single Panel Test Type is manufactured using environmentally friendly materials, where possible. Evans recommend that the equipment is dismantled as far as possible for recycling purposes. Electrical parts and cables should be disposed of as electric waste, in accordance with local legislation.

Serious Incident Statement:

Please immediately report any serious incident involving this product back to Evans (Instruments) Limited, together with the competent authority within your country (MHRA in the UK).

Guarantee

If this product becomes defective due to faulty materials and workmanship, within 12 months from the date of purchase, we guarantee to either replace all defective parts (with the exception of LEDs) or, at our discretion, replace the unit free of charge, provided that:

- The product is returned to us or our authorised repairers with evidence of date of purchase;
- The product has not been subjected to misuse or neglect;
- The product has not sustained any damage through foreign objects, substances or accidents;
- Repairs have not been attempted by anyone other than our own service staff or authorised repair distributors.

This guarantee is offered as an extra benefit and is additional to the customer's statutory rights.

Evans (Instruments) Limited
35 Howlett Way, Thetford, Norfolk, IP24 1HZ
Tel: 01842 766004
e: sales@evansinstruments.co.uk
www: evansinstruments.co.uk

EC REP

VeriSafe Compliance Limited
77 Lower Camden Street, Dublin, Ireland, DO2 XE80

Single Panel Test Type
900(6-L) / 900(6-L/WB) / 900(6-L/H)



Product Information Booklet

(pictures shown for illustration purpose only)

Please read before use and retain for reference
(Original Instructions)

Issue Date : 07/10/2023

Specifications

GS1-GTIN:	5060660330102	
Model:	900(6-L) & 900(6-L/H)	900(6-L/WB)

Dimensions:Unpacked/Packed

Height:	525/640mm	525/640mm
Width:	210/270mm	210/270mm
Depth:	52/110mm	180/220mm
Weight:	1.5/2.3k	2.3/3.4k

Electrical :
(power supply unit) EN60601-1: approved
Input: 100-240v/AC, 50-60Hz, 160-80mA
Output (to Evans unit): 12v/DC, 600mA, 7.2w
LED Strip (12v)

Power Input Lead Length: 2 metres (entry on right-hand side)
Standard Tests: Test Chart (direct, reverse or dual)
Spot

Spares/Parts List

Part Description	***Part No:***
Neo 6, 12v Medical Plug	EDV1960970EVANS

Cleaning Instructions

Disconnect from mains.
This unit may be cleaned by wiping with a damp cloth and a weak solution of detergent/disinfectant. Please note, alcohol-based cleaning fluids should not be used.

Terms of Use

For indoor use only—protect from moisture.
This product is intended for use within a professional optical practice or clinic environment.
We cannot be held responsible for accidents that may occur after installation of equipment.

Instructions for Use

After fitting the test type in the required position:
The power adaptor supplied with the unit must be connected to an adequate wall power point.
The input plug is inserted into the socket in the bottom right-hand side of the test type.
The panel/spot switch can now be operated, to select the required illuminated section.
If a different chart panel is required, turn off the mains power and slide the panel up to remove from the test type. Reverse the above procedure to fit replacement panel.
Replacement illumination strips should not be required for 30-50,000 hours.
The unit should be returned to Evans (Instruments) Limited for replacement LED strips to be fitted.

Symbols

MD Medical Device.

Class II Equipment.

Operating Instructions.

CE The CE mark on this product indicates it has been tested and conforms with the provisions noted within the Medical Device Regulations (EU)2017/745.

UK CA The UKCA mark on this product indicates it has been tested and conforms with the provisions noted within the UK Medical Device Regulations 2002 (as amended).

The symbol adopted by the European Council to represent waste electrical and electronic equipment comprises a crossed-out wheelie bin with or without a single black line underneath the symbol. The black line indicates that goods have been placed on the market after 2005, when the Directive came into force.

Caution.

SN Serial number of Device.

REF Catalogue Number.

Date of Manufacture.

Manufacturer Details.

EC REP European Authorised Representative.