

OMRON C28P NEBULIZER

• **28098 OMRON C28P COMPRESSOR NEBULIZER - NE-C105-E**
C28P nebulizer with Advanced Valve Technology has a breath-enhanced valve system that ensures maximum drug flow during inhalation and no waste of medication during exhalation.

OMRON

It can be used to treat chronic conditions or as relief from cough and cold, with medications or saline solution. Designed for frequent use by family members and medical facilities.

- efficient treatment of lower airway conditions
- long-life compressor
- durable and easy to clean
- **28099 CHILD MASK - spare**
- **28100 ADULT MASK - spare**



28098



TECHNICAL SPECIFICATIONS

Capacity:	min. 2 ml - max. 10 ml
Operating mode:	continuous use
Particle size:	approx. 2.65 µm
Nebulizing rate:	approx 0.12 ml/min
Noise level:	55 dB
Voltage:	230 V - 50 Hz, 140 VA
Size - weight:	180x180x104 mm - 1.3 kg
Safety:	Type BF (applied part), Class II, IP21

STANDARD ACCESSORIES

Nebulizer kit
Child and adult mask
PVC air tube (150 cm)
Mouthpiece set
5 spare air filters
Power cable and adapter
Carrying bag
Multilingual manual GB, FR, IT, ES, DE, NL, TR, RU, Arabic.

OMRON MICROAIR

• **28080 OMRON MICROAIR U100 MESH NEBULIZER - NE-U100-E**



MicroAIR U100 is a small, portable mesh nebulizer that helps to efficiently treat many respiratory conditions of the lower airways, from a simple cough to asthma.

- **portable**, pocket-sized and battery operated, it can be used anywhere, anytime
 - **soundless operation** enables discreet use
 - **360° operation**: the unique design allows to use it in any position, even when lying down in bed, making it more comfortable for children
 - simple **one-button operation**
 - unique **mesh technology** makes it more efficient and more comfortable to use than traditional nebulizers
- Supplied with unit cover, medication container, mesh cap, mask adapter, mouthpiece, adult mask, child mask, soft pouch, 2xAA batteries, quick guide and multilingual manual GB, FR, IT, ES, DE, NL, TR, RU, Arabic.
- **28081 CHILD MASK - spare**
 - **28082 ADULT MASK - spare**

28080

OMRON



TECHNICAL SPECIFICATIONS

Operating mode:	continuous use
Nebulizing rate:	> 0.25 ml/min
Capacity:	1-10 ml
Particle size:	4.5 µm
Noise level:	20 dB
Power supply:	2xAA batteries, 3 VDC
Battery life:	4 hours
Power consumption:	approx 1.2 W
Size - weight:	38x60x130 mm - 120 g
Safety:	type BF (applied part), IP55



28080

INFRARED LAMP - 150 W

• **28644 INFRARED LAMP WITH PHILIPS BULB 150W - stand**

The infrared lamp can provide localized heating through infrared radiations.

The lamp has 360° adjustable irradiation head.

Adjustable height from 60 to 150 cm.

Supplied with Philips bulb.


Size: 31.5 x h 60 to 150 cm. Weight: 6 kg.

Voltage: AC 210-230 V

• **28647 PHILIPS BULB 150W - spare**

28644

Multilingual manual:
GB, FR, IT, ES, PT,
DE, SE, PL, RO, GR.



AIRO MESH NEBULIZERS

• **28077 AIRO MESH NEBULIZER - GB, ES, PT, NL, SE, LT, PL, CZ, HU, HR, RO, GR**
• **28078 AIRO MESH NEBULIZER - IT, FR, DE, DK, NO, FI, LV, EE, SK, SI, BG, Arabic.**

This ultra-compact and lightweight mesh nebuliser is designed for both clinical and home environments.

By converting liquid medication into a fine, gentle mist, it ensures that particles are easy to inhale and absorb, providing quick and effective relief.

- one-button operation with single short press allows users to quickly switch between working modes
- high and adjustable nebulisation rate with two working modes
- low-noise
- automatically turns off when the liquid reservoir is empty
- built-in rechargeable battery
- power indication

Supplied with user manual, type C power cable, mouthpiece, adult and child masks.

TECHNICAL SPECIFICATIONS

Nebulization rate:	level I: ≥ 0.25 ml/min level II: ≥ 0.15 ml/min
Working frequency:	110 KHz, +/- 10%
Particle size:	approx. 5 µm
Operating mode:	continuous
Power supply:	3.7 V lithium battery or power adapter DC 5V
Size:	59.4x48.7x113.8 mm
Protection:	Class II, type BF applied part



28077-28078

HOT-COLD THERMY GEL

• **28593 HOT-COLD THERMY GEL 14x18 cm - box of 25**

• **28594 HOT-COLD THERMY GEL 11x26 cm - box of 25**

Non-toxic and reusable hot/cold gel suitable for medical purposes. These gel packs are microwave and boiling water-safe, making them ideal for heat therapy: lumbago, stiff-neck, rheumatism, hypothermia and whenever hot application is required. When placed in the freezer, they are perfect for cold therapy: haematomas, sprains, injuries, hyperthermia and whenever cold application is required.

Thermy gel 28594 can be used with a cover (28596) to prevent burns. Class I, Medical Device.

• **28596 COVER WITH FASTENING VELCRO 25.5x12.5 cm for 28594**



28593

Multilingual instructions in:

GB*, FR*, IT*, ES*, PT*, DE*, PL, NL, CZ, DK, SE, FI, EE, LT, LV, SK, SI, RO, HU, HR, GR, BG, Arabic. *Also printed on product



28594



28594+28596