Shower and commode chairs

REBOTEC®



www.rebotec.de

Köln / Kiel

Commode chairs

Equipment:

- ergonomically shaped backrest
- swinging armrests enable easy side transfer
- toiletseat with seat cushion
- bucket device removable
- incl. bucket
- nonsliding ferrules for safe standing
- approved up to 130 kg body weight

Köln

ref. no. 340. .. .00

Kiel

- seat height adjustable

ref. no. 347. .. .00

Accessories:

- armrest cushion ref. no. 440.10.32
- holder for disposable bags with lid ref.no. 304.05.10
- disposable bags with lid ref.no. 301.00.20



ref. no. 340.44.00



ref. no. 347.03.00

Type:	Köln	Kiel
depth:	63 cm	63 cm
width:	57 cm	57 cm
seat depth:	45 cm	45 cm
seat height:	52 cm	51 cm - 56 cm
seat width:	46 cm	46 cm
weight:	11,2 kg	11,2 kg
max. body weight:	130 kg	130 kg
ref. no.	340. .. .00	347. .. .00

Colour	Colour code
grey/grey	03
grey/rubin-red	44
grey/blue	54

Technical changes and changes in colour reserved.

protected by patent

CE

Bonn

Commode chair



ref. no. 343.54.00

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- toiletseat with seat cushion
- bucket device removable
- incl. bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake
- approved up to 130 kg body weight

ref. no. 343. .. .00

Accessories:

- armrest cushion
  ref. no. 440.10.32
- safety belt
  ref. no. 322.30.20
- IV pole
  ref. no. 303.62
- holder for disposable bags with lid
  ref.no. 304.05.10
- disposable bags with lid
  ref.no. 301.00.20



ref. no. 343.10.22
Also available with hygiene opening!

Colour	Colour code
grey/grey	03
grey/black	10
terracotta/rubin-red	33
grey/rubin-red	44
grey/blue	54
colourful	97

Type:	Bonn
depth:	87 cm
width:	56 cm
seat depth:	45 cm
seat height:	53 cm
seat width:	46 cm
weight:	14,5 kg
toilet-height with bucket device:	44 cm
toilet-height without bucket device:	47 cm
max. body weight:	130 kg
ref. no.	343. .. .00

Technical changes and changes in colour reserved.

protected by patent

Hamburg

Commode chair, height adjustable

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- toiletseat with seat cushion
- bucket device removable
- incl. bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake
- approved up to 130 kg body weight

ref. no. 356. .. .00

Accessories:

- armrest cushion ref. no. 440.10.32
- safety belt ref. no. 322.30.20
- IV pole ref. no. 303.62
- holder for disposable bags with lid ref.no. 304.05.10
- disposable bags with lid ref.no. 301.00.20



ref. no. 356.44.00

Type:	Hamburg
depth:	87 cm
width:	56 cm
seat depth:	45 cm
seat height:	54 cm - 59 cm
seat width:	46 cm
toilet-height with bucket device:	45 cm - 50 cm
toilet-height without bucket device:	48 cm - 53 cm
weight:	17 kg
max. body weight:	130 kg
ref. no.	356. .. .00

Colour	Colour code
grey/grey	03
grey/rubin-red	44
grey/blue	54



Technical changes and changes in colour reserved.

protected by patent

CE

Boston

Commode chair up to 150 kg body weight



ref. no. 342.03.00

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- toiletseat with seat cushion
- bucket device removable
- incl. bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake
- approved up to 150 kg body weight

ref. no. 342. .. .00

Accessories:

- armrest cushion
  ref. no. 440.10.32
- safety belt
  ref. no. 322.30.20
- IV pole
  ref. no. 303.62
- holder for disposable bags with lid
  ref.no. 304.05.10
- disposable bags with lid
  ref.no. 301.00.20

Up to 150 kg!

Colour	Colour code
grey/grey	03
grey/blue	54

Type:	Boston
depth:	89 cm
width:	62 cm
seat depth:	45 cm
seat height:	53 cm
seat width:	52 cm
toilet-height with bucket device:	44 cm
toilet-height without bucket device:	47 cm
weight:	16,2 kg
max. body weight:	150 kg
ref. no.	342. .. .00

Technical changes and changes in colour reserved.

protected by patent

CE

Zwickau / Würzburg

Shower chairs

Equipment:

- ergonomically shaped backrest
- swinging armrests enable easy side transfer
- nonholstered seat with hygiene opening
- nonsliding ferrules for safe standing
- approved up to 130 kg body weight

Zwickau

ref. no. 332. .. .20

Würzburg

- seat height adjustable

ref. no. 333. .. .20

Accessories:

- bucket
  ref. no. 301. ..
- bucket device
  ref. no. 309.00.00
- seat cushion
  ref.no. 440. .. .62
- armrest cushion
  ref. no. 440.10.32
- soft seat with hygiene opening
  ref.no. 440.10.04



ref. no. 332.54.20

New colours!



ref. no. 333.54.20

Type:	Zwickau	Würzburg
depth:	63 cm	63 cm
width:	57 cm	57 cm
seat depth:	45 cm	45 cm
seat height:	51 cm	50 cm - 55 cm
seat width:	46 cm	46 cm
weight:	9,2 kg	9,2 kg
max. body weight:	130 kg	130 kg
ref. no.	332. .. .20	333. .. .20

Colour	Colour code
grey/grey	03
grey/black	10
grey/rubin-red	44
grey/blue	54

Technical changes and changes in colour reserved.

protected by patent

CE

Erfurt

Shower chair with 5-inch castors



ref. no. 334.10.20

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- nonholstered seat with hygiene opening
- height adjustable, swinging and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake
- approved up to 130 kg body weight

Erfurt

ref. no. 334. .. .20

Erfurt

- incl. bucket device and bucket

ref. no. 334. .. .21

Erfurt

- incl. bucket device, bucket and seat cushion

ref. no. 334. .. .22

Accessories:

- armrest cushion ref. no. 440.10.32
- safety belt ref. no. 322.30.20
- IV pole ref. no. 303.62



ref. no. 334.44.22

Colour	Colour code
grey/grey	03
grey/black	10
grey/rubin-red	44
grey/blue	54

Type:	Erfurt
depth:	87 cm
width:	56 cm
seat depth:	45 cm
seat height:	53 cm
seat width:	46 cm
toilet-height with bucket device:	43 cm
toilet-height without bucket device:	47 cm
weight:	11,9 kg
max. body weight:	130 kg
ref. no.:	334. .. . ..

Technical changes and changes in colour reserved.

protected by patent

CE

Augsburg

Shower-commode chair for children

Equipment:

- aluminium back with PUR backrest
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening
- incl. bucket
- height adjustable swinging and detachable foot support
- 4 swivelling castors with parking/direction brake
- approved up to 80 kg body weight

ref. no. 339. .. .97

Types:

3 inch ref. no. 339.03.97
4 inch ref. no. 339.04.97
24 inch ref. no. 339.24.97

Accessories:

- height adjustable headrest ref. no. 444.03.39
- safety belt ref. no. 322.30.10
- seat width and depth adjustment ref. no. 444.00.97
- armrest cushion ref. no. 440.10.32



ref. no. 339.04.97

Also available as 24-inch selfpropelling chair!



ref. no. 339.24.97

Type:	3 inch	4 inch	24 inch
depth:	87 cm	87 cm	97 cm
width:	53 cm	53 cm	66 cm
seat depth:	35 cm	35 cm	35 cm
seat height:	51 cm	54 cm	55 cm
seat width:	40 cm	40 cm	40 cm
weight:	17,1 kg	17,3 kg	24,8 kg
toilet-height with bucket device:	40 cm	43 cm	44 cm
toilet width:	33 cm	33 cm	32 cm
max. body weight:	80 kg	80 kg	80 kg
ref. no.	339.03.97	339.04.97	339.24.97

Technical changes and changes in colour reserved.

protected by patent

CE

Seat width and depth adjustment

for shower-commode chair Augsburg



ref. no. 339.04.97
shown with height adjustable head-rest
and seat width and depth adjustment


Equipment:

- backrest adjustment
- adjustable armrests
- safety belt

ref. no. 444.00.97

Accessories:

- height adjustable headrest
  ref. no. 444.03.39





Up to 80 kg!

Type:	Seat width and depth adjustment
seat width:	22 cm - 32 cm
seat depth:	25 cm - 29 cm
weight:	3,9 kg
ref. no.	444.00.97

Technical changes and changes in colour reserved.

protected by patent

CE

Konstanz / Ulm

Shower-commode chairs

Equipment:

- ergonomically shaped backrest
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening (available in rear or in front)
- incl. bucket device and bucket
- nonsliding ferrules for safe standing
- approved up to 130 kg body weight

Konstanz

ref. no. 341. .. .00

Ulm

- seat height adjustable

ref. no. 348. .. .00

Accessories:

- armrest cushion ref. no. 440.10.32
- shower seat without hygiene opening ref. no. 440.10.24



ref. no. 341.54.00



ref. no. 348.30.00

Type:	Konstanz	Ulm
depth:	63 cm	63 cm
width:	57 cm	57 cm
seat depth:	45 cm	45 cm
seat height:	54 cm	52 cm - 57 cm
seat width:	46 cm	46 cm
weight:	11,8 kg	11,8 kg
max. body weight:	130 kg	130 kg
ref. no.	341. .. .00	348. .. .00

Colour	Colour code
grey/grey	03
grey/fire-red	30
grey/rubin-red	44
grey/blue	54

Technical changes and changes in colour reserved.

protected by patent

CE

Berlin

Shower-commode chair with 5-inch castors



ref. no. 344.33.00

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening (available in rear or in front)
- incl. bucket device and bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake
- approved up to 130 kg body weight

ref. no. 344. .. .00

Accessories:

- IV pole
  ref. no. 303.62
- safety belt
  ref. no. 322.30.20
- armrest cushion
  ref. no. 440.10.32
- shower seat without hygiene opening
  ref. no. 440.10.24

Colour	Colour code
grey/grey	03
grey/fire-red	30
terracotta/rubin-red	33
grey/rubin-red	44
grey/blue	54
colourful	97

Type:	Berlin
depth:	87 cm
width:	56 cm
seat depth:	45 cm
seat height:	55 cm
seat width:	46 cm
weight:	15,9 kg
toilet-height with bucket device:	43 cm
toilet-height without bucket device:	47 cm
max. body weight:	130 kg
ref. no.	344. .. .00

Technical changes and changes in colour reserved.

protected by patent

CE

Frankfurt

Shower-commode chair, height adjustable

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening (available in rear or in front)
- incl. bucket device and bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake
- approved up to 130 kg body weight

ref. no. 355. .. .00

Accessories:

- IV pole
  ref. no. 303.62
- safety belt
  ref. no. 322.30.20
- armrest cushion
  ref. no. 440.10.32
- shower seat without hygiene opening
  ref. no. 440.10.24



ref. no. 355.54.00

Type:	Frankfurt
depth:	87 cm
width:	56 cm
seat depth:	45 cm
seat height:	56 cm - 61 cm
seat width:	46 cm
toilet-height with bucket device:	44 cm - 49 cm
toilet-height without bucket device:	48 cm - 53 cm
weight:	18 kg
max. body weight:	130 kg
ref. no.	355. .. .00

Colour	Colour code
grey/grey	03
grey/fire-red	30
grey/rubin-red	44
grey/blue	54



Technical changes and changes in colour reserved.

protected by patent

CE

New York

Shower-commode chair up to 150 kg



ref. no. 369.54.00

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening (available in rear or in front)
- incl. bucket device and bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake
- approved up to 150 kg body weight

ref. no. 369. .. .00

Accessories:

- IV pole
  ref. no. 303.62
- safety belt
  ref. no. 322.30.20
- armrest cushion
  ref. no. 440.10.32
- shower seat without hygiene opening
  ref. no. 440.10.24

Up to 150 kg!

Colour	Colour code
grey/fire-red	30
grey/blue	54

Type:	New York
depth:	89 cm
width:	62 cm
seat depth:	45 cm
seat height:	55 cm
seat width:	52 cm
toilet-height with bucket device:	43 cm
toilet-height without bucket device:	47 cm
weight:	16,6 kg
max. body weight:	150 kg
ref. no.	369. .. .00

Technical changes and changes in colour reserved.

protected by patent

CE

Denver / Atlanta

Shower-commode chairs up to 175 kg

Equipment:

- aluminium back with PUR backrest
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening
- incl. bucket device and bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors, in the rear with parking/direction brake
- approved up to 175 kg body weight

Denver
ref. no. 360. .. .00

Atlanta
ref. no. 361. .. .00

Accessories:

- armrest cushion ref. no. 440.10.32
- safety belt ref.no. 322.30.20



ref. no. 360.54.00



ref. no. 361.30.00

Up to 175kg!

Type:	Denver	Atlanta
depth:	90 cm	90 cm
width:	68 cm	74 cm
seat depth:	48 cm	48 cm
seat height:	57 cm	57 cm
seat width:	58 cm	64 cm
toilet-height with bucket device:	45 cm	45 cm
toilet-height without bucket device:	48 cm	48 cm
weight:	19 kg	20 kg
max. body weight:	175 kg	175 kg
ref. no.	360. .. .00	361. .. .00

Colour	Colour code
grey/fire-red	30
grey/blue	54

Technical changes and changes in colour reserved.

protected by patent

CE

Chicago / Dallas

Shower-commode chairs up to 200 kg



ref. no. 358.30.00



ref. no. 359.54.00

Equipment:

- aluminium back with PUR backrest
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening
- incl. bucket device and bucket
- height adjustable, swinging and detachable foot support
- 4 swivelling castors, in the rear with parking/direction brake
- approved up to 200 kg body weight

Chicago
ref. no. 358. .. .00

Dallas
ref. no. 359. .. .00

Accessories:

- armrest cushion ref. no. 440.10.32
- safety belt ref. no. 322.30.20

Up to 200 kg!

Colour	Colour code
grey/fire-red	30
grey/blue	54

Type:	Chicago	Dallas
depth:	90 cm	90 cm
width:	68 cm	74 cm
seat depth:	48 cm	48 cm
seat height:	56 cm	56 cm
seat width:	58 cm	64 cm
toilet-height with bucket device:	44 cm	44 cm
toilet-height without bucket device:	47 cm	47 cm
weight:	21 kg	22 kg
max. body weight:	200 kg	200 kg
ref. no.	358. .. .00	359. .. .00

Technical changes and changes in colour reserved.

protected by patent

CE

Genf

Selfpropelling shower-commode chair

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- height adjustable, swinging and detachable foot support
- 24-inch driving wheels with grip wheel
- approved up to 130 kg body weight

Genf

- upholstered seat with hygiene opening
- incl. bucket device and bucket

ref. no. 346. .. .10

Genf

- seat with hygiene opening
- incl. bucket device and bucket

ref. no. 346. .. .20

Genf

- toiletseat
- incl. bucket device and bucket

ref. no. 346. .. .30

Genf

- upholstered seat without hygiene opening

ref.no. 346. .. .40



ref. no. 346.30.10

Up to 130 kg!

Colour	Colour code
grey/grey	03
grey/black	10
grey/fire-red	30
grey/rubin-red	44
grey/blue	54

Type:	Genf	Genf	Genf
depth:	95 cm	95 cm	95 cm
width:	71 cm	71 cm	71 cm
seat depth:	45 cm	45 cm	45 cm
seat height:	55 cm	52 cm	52 cm
seat width:	46 cm	46 cm	46 cm
toilet-height with bucket device:	43 cm	43 cm	43 cm
max. toilet width:	38 cm	38 cm	38 cm
weight:	24 kg	22 kg	23 kg
max. body weight:	130 kg	130 kg	130 kg
ref. no.	346. .. .10	346. .. .20	346. .. .30

Technical changes and changes in colour reserved.

protected by patent

CE

Miami

Selfpropelling shower-commode chair



ref. no. 345.10.10

Equipment:

- ergonomically shaped backrest with integrated grip
- swinging armrests enable easy side transfer
- height adjustable, swinging and detachable foot support
- 24-inch driving wheels with grip wheel
- approved up to 150 kg body weight

Miami

- upholstered seat with hygiene opening
- incl. bucket device and bucket

ref. no. 345. .. .10

Miami

- seat with hygiene opening
- incl. bucket device and bucket

ref. no. 345. .. .20

Miami

- toiletseat
- incl. bucket device and bucket

ref. no. 345. .. .30

Miami

- upholstered seat without hygiene opening

ref.no. 345. .. .40

Up to 150 kg!

Type:	Miami	Miami	Miami
depth:	95 cm	95 cm	95 cm
width:	77 cm	77 cm	77 cm
seat depth:	45 cm	45 cm	45 cm
seat height:	55 cm	52 cm	52 cm
seat width:	52 cm	52 cm	52 cm
toilet-height with bucket device:	43 cm	43 cm	43 cm
weight:	27 kg	24 kg	25 kg
max. body weight:	150 kg	150 kg	150 kg
ref. no.	345. .. .10	345. .. .20	345. .. .30

Colour	Colour code
grey/grey	03
grey/black	10
grey/fire-red	30
grey/rubin-red	44
grey/blue	54

Technical changes and changes in colour reserved.

protected by patent

CE

Dubai

Toilet chair with washing seat

- shower-toilet-seat for optimal cleaning
- swinging armrests enable easy side transfer
- 4 swivelling castors 5-inch with parking/direction brake

Equipment:

- Powerwash-Function: bubbling jet increases the efficiency of cleaning
- massage by a pulsating water jet promotes blood circulation
- washing position, water pressure and temperature selection
- adjustable hot air dryer
- heated seat
- programming of individual user profiles
- Enema-Function

ref. no. 329.05. ..



ref. no. 329.05.30



Up to 150 kg!

Odour extraction and warm-water -cleaning!





Colour	Colour code
grey/grey	03
grey/fire-red	30
grey/blue	54

Type:	Dubai
depth:	92 cm
width:	63 cm
seat height:	56 cm
seat width:	52 cm
weight:	25 kg
max. body weight:	150 kg
ref. no.	329.05. ..

Technical changes and changes in colour reserved.

protected by patent

CE

Malta

Shower and tub-seat chair



ref. no. 367.54.30

Up to 130 kg!

- The tub-seat chair for seated bathing, washing and hygiene of the lower abdominal area.

Equipment:

- aluminium back with PUR backrest
- swinging armrests enable easy side transfer
- upholstered seat with hygiene opening
- height adjustable, swinging and detachable foot support
- 4 swivel castors 3-inch with parking/direction brake
- approved up to 130 kg body weight

Tub dimensions:

- tub height: 17 cm
- tub length: 56 cm
- tub width: 42 cm
- volume: 6-15 l

ref. no. 367.54.30 grey/blue

Accessories:

- armrest cushion ref. no. 440.10.32





Colour	Colour code
grey/blue	54

Type:	Malta
depth:	90 cm
width:	74 cm
seat height:	48 cm
seat width:	64 cm
seat depth:	48 cm
tub rim height:	42 cm
weight:	19,3 kg
max. body weight:	130 kg
ref. no.	347. .. .30







Technical changes and changes in colour reserved.

protected by patent

CE

Oslo

Transport and shower chair

Equipment:

- aluminium back with PUR backrest
- swinging armrests enable easy side transfer
- upholstered seat
- height adjustable, swinging and detachable foot support
- 4 swivelling castors, in the rear with parking/direction brake
- approved up to 150 kg body weight

ref. no. 362. .. .00

Accessories:

- IV pole
  ref. no. 303.62
- safety belt
  ref. no. 322.30.20
- armrest cushion
  ref. no. 440.10.32



ref. no. 362.54.00

Convenient handle and comfortable easy-to-clean seat!



Up to 150 kg!

Colour	Colour code
grey/fire-red	30
grey/blue	54

Type:	Oslo
depth:	106 cm
width:	64 cm
seat height:	56 cm
seat width:	52 cm
seat depth:	45 cm
weight:	16,5 kg
max. body weight:	150 kg
ref. no.	362. .. .00

Technical changes and changes in colour reserved.

protected by patent

CE

Accessories



Colour	Colour code
grey	00
rubin-red	44
brown	48

1. Bucket

- suitable for all commode chairs and shower chairs

ref. no. 301. ..



2. Holder for disposable bags with lid

- incl. 3 disposable bags with lid

ref.no. 304.05.10



3. Disposable bags with lid

ref.no. 301.00.20



ref. no. 307.00

Colour	Colour code
grey	00
rubin-red	44
blue	54

4. Toilet seat

- consisting of two parts, bucket removeable upwardly and rearwardly

ref.no. 307. ..

Technical changes and changes in colour reserved.

protected by patent

CE

Accessories

5. Shower seat without hygiene opening

- suitable for all models
  with seat width 46 and 52 cm

ref. no. 440.10.24



6. Splash-protection

- suitable for soft seats with hygiene opening
  and seat width 46 and 52 cm

ref.no. 440.10.40



7. Safety belt

size S

ref. no. 322.30.10

size M/L

ref. no. 322.30.20



8. Cushion for safety belt

ref. no. 322.10.50




Technical changes and changes in colour reserved.

protected by patent

CE

Accessories



9. Armrest cushion

- suitable for all models

ref. no. 440.10.32

10. IV pole

- suitable for all models (except 24-inch, bariatric models and Phoenix models)

ref. no. 440.00.93

Technical changes and changes in colour reserved.

protected by patent

CE

www.rebotec.de

Distributed by:

06-2022

11702-19306

REBOTEC® Rehabilitationsmittel GmbH

D-49610 Quakenbrück · Artlandstraße 57 – 59 · Fon +49 (0) 5431 9416-0 · Fax +49 (0) 5431 9416-66

www.rebotec.de · info@rebotec.de