Bristol Maid™
Hospital Metalcraft™ Limited

Email:
sales@bristolmaid.com
Phone: +44 (0)1258
451338
Fax: +44 (0)1258 455056

Pharmacy Trolleys - Trays

Product Images











Description

FEATURES

- Trolleys suitable for the storage & dispensing of drugs & medicines using individual patient trays
- Two locking options...
  - High Security Bolt Lock - Supplied with two keys
  - Electronic Push Button - Programmable four-digit code
- Trays supported on adjustable shelves, sizes (w x d x h)...
  - LA - 86 x 268 x 57mm (aluminium)
  - A - 85 x 181 x 54mm (plastic)
  - C - 129 x 181 x 54mm (plastic)
  - LP - 115 x 266 x 95mm (plastic)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- Adjustable door shelves...

    - Single Door - Three door shelves
    - Double Door - Six door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building - see optional MT/1/6F

STANDARD COLOURS

- Paint - Grey White
- Work Surface - Spearmint

Product Options

Code	Description	Overall Dimensions (W x D x H)
PTS/HSB/LA16	Sixteen LA Trays - High Security Bolt Lock	511 x 542 x 1072mm
PTD/HSB/A18C12	Eighteen A & Twelve C Trays - High Security Bolt Lock	711 x 568 x 1072mm
PTD/HSB/LP24	Twenty-Four LP Trays - High Security Bolt Lock	711 x 568 x 1072mm
PTS/EPB/LA16	Sixteen LA Trays - Electronic Push Button Lock	511 x 542 x 1072mm
PTD/EPB/A18C12	Eighteen A & Twelve C Trays - Electronic Push Button Lock	711 x 568 x 1072mm
PTD/EPB/LP24	Twenty-Four LP Trays - Electronic Push Button Lock	711 x 568 x 1072mm