EnSafe.One®
A world open to all


The **HULK Bariatric** evacuation mattress allows the safe evacuation of a bariatric person in emergency situations.

evacuation is carried out by two people

designed for transport down stairs and on flat surfaces

an ideall solution for bariatric patients who are bedridden or post-operative

- **Model:** HULK Bariatric
- **Product Class:** Class I non-invasive medical device
- **Dimensions:** 205 x 85 x 7 cm [H x W x D]
- **Material:** flame-retardant rubberized polyester, waterproof and water-resistant, polyamide / foam / velour / PVC. Includes 1 foam insert
- **Securing straps with buckle:** 3 pcs.
- **Securing flaps:** 3 pcs.
- **Transport handles:** 2 pcs (each with 3 different lengths)
- **Side handles with reflective elements:** 12 pcs.
- **Weight:** 7 kg (without cover)
- **Maximum load capacity:** up to 300 kg
- **Additional elements:** detachable, adjustable head support with strap; Tuplex rigid plate under the feet
- **Storage:** carry bag or mountable cover

















www.ensafe.one

HULK BARIATRIC

Manufacturer: EnSafe.One A world open to all

This is a medical device.
Use it according to the instructions or label.

ver.: 05_05_2026