
KERN **MBC-A08**

Rechargeable battery 6×1.2 V 2000 mA 7,2 V (NiMH, 6 x 1.2 V 2000 mA)



Category	
Brand	KERN
Product category	Balance
Product group	Rechargeable battery

Construction	
Dimension (W×D×H)	30×150×15 mm

Power Supply	
Battery / accumulator type	NiMH
Battery connection	Plug-in terminals
Battery capacity	2.000 mAh
Battery voltage	7,2 V
Battery capacity (Wh)	14,4 Wh
Rechargeable battery	6×1.2 V 2000 mA

Packing & Shipping	
Delivery time	1 d
Dimensions packaging (W×D×H)	35×155×20 mm
Shipping method	Parcel service
Net weight approx.	0,15 kg
Gross weight approx.	0,20 kg
Shipping weight	0,163 kg

Product information	
GTIN/EAN number	4045761177635
Trade number	85075000

Pictograms

STANDARD

